


## Medtronic Shunt Valves

<http://www.medtronic.com/neurosurgery/shunts.html>

Medtronic Strata II Valve (a programmable valve, different settings are depicted; “P/L” stands for Performance Level; see Medtronic website for pressure/flow information)



**P/L 0.5**


**P/L 1.0**

**P/L 1.5**


**P/L 2.0**

**P/L 2.5**


Medtronic Strata II Valve  
Image depicts valve set at P/L 1.5




Medtronic Strata Valve, Small  
Image depicts valve set at P/L 1.0



Medtronic Delta<sup>®</sup> Valve, Performance Level 1.0



Medtronic Delta® Burr Hole Valve, Performance Level 2.0



Medtronic Delta® Valve Performance Level Code  
Arrows depict direction of flow



Level 0.5



Level 1.0



Level 1.5




Level 2.0

Medtronic CSF-Flow Control valves  
No xray images available. A photograph is shown.



Medtronic CSF- Flow Control Valve, Pressure Level Code  
Arrows depict direction of flow

High Pressure:



Medium Pressure:



Low Pressure:



Low-Low Pressure:



Strata<sup>®</sup> and Delta<sup>®</sup> are registered trademarks of Medtronic, Inc.